

# Sock Aid

1. Designed for those with limited hand function or back problems, wheelchair users, or those recovering from hip or knee surgery.
2. reduce waist bending when wearing socks.
3. Easy to use.

## User manual:

- 1、 Slide your sock over the sleeve of the sock aid.
- 2、 Holding the cloth strap, lower the sock aid to the floor
- 3、 Slip your foot into the opening of the sock and pull the strap up until the sock fits comfortably on your foot.
- 4、 Pull further on the strap to remove the sock aid.

## Wash instructions:

Hand wash with mild detergent, hang dry.

